

Product Datasheet

Human HS6ST2 ELISA Kit (orb440102)

Description	Human Heparan Sulfate-6-O-Sulfotransferase 2 (HS6ST2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q96MM7
Sensitivity	0.066ng/mL
Entrez	90161
Expiration Date	Please enquire.