

Product Datasheet

Human HSD11b1 ELISA Kit (orb440104)

Description	Human 11-Beta-Hydroxysteroid Dehydrogenase Type 1 (HSD11b1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P28845
Sensitivity	0.57ng/mL
Entrez	3290
Expiration Date	Please enquire.