

## Product Datasheet

### Human HSP90aA1 ELISA Kit (orb440123)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Heat Shock Protein 90kDa Alpha A1 (HSP90aA1) Sandwich ELISA Kit                                   |
| <b>Reactivity</b>        | Human                                                                                                   |
| <b>Range</b>             | 1.56-100ng/mL                                                                                           |
| <b>Form/Appearance</b>   | Traditional                                                                                             |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.                |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| <b>Uniprot ID</b>        | <b>P07900</b>                                                                                           |
| <b>Sensitivity</b>       | 0.55ng/mL                                                                                               |
| <b>Entrez</b>            | <b>3320</b>                                                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                                                         |