

Product Datasheet

Bovine HSP90aB1 ELISA Kit (orb440126)

Description	Bovine Heat Shock Protein 90kDa Alpha B1 (HSP90aB1) Sandwich ELISA Kit
Reactivity	Bovine
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q76LV1
Sensitivity	0.71ng/mL
Entrez	767874
Expiration Date	Please enquire.