

Product Datasheet

Mouse HSP90b1 ELISA Kit (orb440131)

Description	Mouse Heat Shock Protein 90kDa Beta 1 (HSP90b1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P08113
Sensitivity	0.064ng/mL
Entrez	22027
Expiration Date	Please enquire.