

Product Datasheet

Rat HSPA5 ELISA Kit (orb440145)

Description	Rat Heat Shock 70kDa Protein 5 (HSPA5) Sandwich ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P06761
Sensitivity	0.27ng/mL
Entrez	25617
Expiration Date	Please enquire.