

Product Datasheet

Rat HSPA8 ELISA Kit (orb440148)

Description	Rat Heat Shock 70kDa Protein 8 (HSPA8) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P63018
Sensitivity	0.103ng/mL
Entrez	24468
Expiration Date	Please enquire.