

Product Datasheet

Chicken HSPb2 ELISA Kit (orb440153)

Description	Chicken Heat Shock Protein Beta 2 (HSPb2) Sandwich ELISA Kit
Reactivity	Gallus
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	R4GG05
Sensitivity	5.7pg/mL
Entrez	414136
Expiration Date	Please enquire.