

Product Datasheet

Human HSPb3 ELISA Kit (orb440156)

Description	Human Heat Shock Protein Beta 3 (HSPb3) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q12988
Sensitivity	0.29ng/mL
Entrez	8988
Expiration Date	Please enquire.