

Product Datasheet

Mouse HSPb6 ELISA Kit (orb440158)

Description	Mouse Heat Shock Protein Beta 6 (HSPb6) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q5EBG6
Sensitivity	0.067ng/mL
Entrez	243912
Expiration Date	Please enquire.