

Product Datasheet

Rat HSPb6 ELISA Kit (orb440159)

Description	Rat Heat Shock Protein Beta 6 (HSPb6) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P97541
Sensitivity	0.068ng/mL
Entrez	192245
Expiration Date	Please enquire.