

Product Datasheet

Human HSPBP1 ELISA Kit (orb440164)

Description	Human Heat Shock 70kDa Binding Protein 1 (HSPBP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9NZL4
Sensitivity	0.105ng/mL
Entrez	23640
Expiration Date	Please enquire.