

Product Datasheet

Human HSPG2 ELISA Kit (orb440169)

Description	Human Heparan Sulfate Proteoglycan 2 (HSPG2) Sandwich ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P98160
Sensitivity	49pg/mL
Entrez	3339
Expiration Date	Please enquire.