

Product Datasheet

Human HTN1 ELISA Kit (orb440173)

Description	Human Histatin 1 (HTN1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P15515
Sensitivity	0.053ng/mL
Entrez	3346
Expiration Date	Please enquire.