

Product Datasheet

Human HTRA4 ELISA Kit (orb440180)

Description	Human HtrA Serine Peptidase 4 (HTRA4) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P83105
Sensitivity	0.32ng/mL
Entrez	203100
Expiration Date	Please enquire.