

## Product Datasheet

### Human HYOU1 ELISA Kit (orb440187)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Hypoxia Up Regulated 1 (HYOU1)<br>Sandwich ELISA Kit                               |
| <b>Reactivity</b>        | Human                                                                                    |
| <b>Range</b>             | 0.156-10ng/mL                                                                            |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Uniprot ID</b>        | <b>Q9Y4L1</b>                                                                            |
| <b>Sensitivity</b>       | 0.061ng/mL                                                                               |
| <b>Entrez</b>            | <b>10525</b>                                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                                          |