

Product Datasheet

Human ICAM1 ELISA Kit (orb440196)

Description	Human Intercellular Adhesion Molecule 1 (ICAM1) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05362
Sensitivity	31pg/mL
Entrez	3383
Expiration Date	Please enquire.