

Product Datasheet

Human ICAM5 ELISA Kit (orb440205)

Description	Human Intercellular Adhesion Molecule 5 (ICAM5) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UMF0
Sensitivity	0.068ng/mL
Entrez	7087
Expiration Date	Please enquire.