

Product Datasheet

Human ICOS ELISA Kit (orb440207)

Description	Human Inducible T-Cell Co Stimulator (ICOS) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine, cell lysates and other biological fluids.
Sample Types	serum, plasma, urine, cell lysates and other biological fluids.
Uniprot ID	Q9Y6W8
Sensitivity	6.2pg/mL
Entrez	29851
Expiration Date	Please enquire.