

Product Datasheet

Human IDS ELISA Kit (orb440219)

Description	Human Iduronate-2-Sulfatase (IDS) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P22304
Sensitivity	0.057ng/mL
Entrez	3423
Expiration Date	Please enquire.