

Product Datasheet

Mouse INSL3 ELISA Kit (orb440260)

Description	Mouse Insulin Like Protein 3 (INSL3) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	009107
Sensitivity	6.3pg/mL
Entrez	16336
Expiration Date	Please enquire.