

Product Datasheet

Mouse IQGAP2 ELISA Kit (orb440277)

Description	Mouse IQ Motif Containing GTPase Activating Protein 2 (IQGAP2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q3UQ44
Sensitivity	0.062ng/mL
Entrez	544963
Expiration Date	Please enquire.