

Product Datasheet

Human ITIH1 ELISA Kit (orb440312)

Description	Human Inter Alpha-Globulin Inhibitor H1 (ITIH1) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P19827
Sensitivity	15.7ng/mL
Entrez	3697
Expiration Date	Please enquire.