

Product Datasheet

Human ITIH5 ELISA Kit (orb440317)

Description	Human Inter Alpha-Globulin Inhibitor H5 (ITIH5) Sandwich ELISA Kit
Reactivity	Human
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q86UX2
Sensitivity	33ng/mL
Entrez	80760
Expiration Date	Please enquire.