

Product Datasheet

Human JAM3 ELISA Kit (orb440335)

Description	Human Junctional Adhesion Molecule 3 (JAM3) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9BX67
Sensitivity	6.2pg/mL
Entrez	83700
Expiration Date	Please enquire.