

Product Datasheet

CD4 antibody (orb44034)

Description

Mouse monoclonal antibody to CD4

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IP

Immunogen

Human thymocytes and T lymphocytes.

Target

CD4

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Immunoprecipitation: Excellent. Flow cytometry: Recommended dilution: 3 µg/ml. Although it has not been tested rigorously, following data suggest that the antibody MEM-115 is a low-affinity antibody: its binding to T cells increases at elevated temperature, monovalent Fab fragments essentially do not bind to T cells.

Isotype

Mouse IgG2a

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

P01730

Entrez


920

Dilution Range


Immunoprecipitation: Excellent. Flow cytometry: Recommended dilution: 3 µg/ml. Although it has not been tested rigorously, following data suggest that the antibody MEM-115 is a low-affinity antibody: its binding to T cells increases at elevated temperature, monovalent Fab fragments essentially do not bind to T cells.

Expiration Date

12 months from date of receipt.



Flow cytometry analysis (surface staining)



Immunoprecipitation of human CD4 from thymocytes