

Product Datasheet

Human KAZALD1 ELISA Kit (orb440344)

Description	Human Kazal Type Serine Protease Inhibitor Domain Containing Protein 1 (KAZALD1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q96I82
Sensitivity	0.115ng/mL
Entrez	81621
Expiration Date	Please enquire.