

Product Datasheet

Human KCC4 ELISA Kit (orb440346)

Description	Human Potassium Chloride Cotransporters 4 (KCC4) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Y666
Sensitivity	0.061ng/mL
Entrez	10723
Expiration Date	Please enquire.