

Product Datasheet

Human Kim1 ELISA Kit (orb440361)

Description	Human Kidney Injury Molecule 1 (Kim1) Sandwich ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, urine, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, urine, cell culture supernatants and other biological fluids.
Uniprot ID	Q96D42
Sensitivity	11.6pg/mL
Entrez	26762
Expiration Date	Please enquire.