

Product Datasheet

Human Klotho ELISA Kit (orb440380)

Description	Human Klotho (KL) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9UEF7
Sensitivity	0.058ng/mL
Entrez	9365
Expiration Date	Please enquire.