

Product Datasheet

Human KLK14 ELISA Kit (orb440386)

Description	Human Kallikrein 14 (KLK14) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9P0G3
Sensitivity	0.044ng/mL
Entrez	43847
Expiration Date	Please enquire.