

Product Datasheet

Human KLK6 ELISA Kit (orb440399)

Description	Human Kallikrein 6 (KLK6) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cerebrospinal fluid, breast milk and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cerebrospinal fluid, breast milk and other biological fluids.
Uniprot ID	Q92876
Sensitivity	0.126ng/mL
Entrez	5653
Expiration Date	Please enquire.