

Product Datasheet

Mouse KLK6 ELISA Kit (orb440400)

Description	Mouse Kallikrein 6 (KLK6) Sandwich ELISA Kit
Reactivity	Mouse
Range	1.56-100pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, cerebrospinal fluid, breast milk and other biological fluids.
Sample Types	serum, plasma, cerebrospinal fluid, breast milk and other biological fluids.
Uniprot ID	Q91Y82
Sensitivity	0.63pg/mL
Entrez	19144
Expiration Date	Please enquire.