

Product Datasheet

Rat KNG1 ELISA Kit (orb440413)

Description	Rat Kininogen 1 (KNG1) Sandwich ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q5PQU1
Sensitivity	11.8pg/mL
Entrez	24903
Expiration Date	Please enquire.