

Product Datasheet

Human L1CAM ELISA Kit (orb440443)

Description	Human L1-Cell Adhesion Molecule (L1CAM) Sandwich ELISA Kit
Reactivity	Human
Range	93.75-6000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P32004
Sensitivity	32.15pg/mL
Entrez	3897
Expiration Date	Please enquire.