

## Product Datasheet

### Human LAIR2 ELISA Kit (orb440449)

|                          |                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Leukocyte Associated<br>Immunoglobulin Like Receptor 2 (LAIR2)<br>Sandwich ELISA Kit  |
| <b>Reactivity</b>        | Human                                                                                       |
| <b>Range</b>             | 15.625-1000pg/mL                                                                            |
| <b>Form/Appearance</b>   | Traditional                                                                                 |
| <b>Storage</b>           | All reagents should be kept according to the<br>instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                       |
| <b>Application notes</b> | serum, plasma, urine, tissue homogenates<br>and other biological fluids.                    |
| <b>Sample Types</b>      | serum, plasma, urine, tissue homogenates<br>and other biological fluids.                    |
| <b>Uniprot ID</b>        | <b>Q6ISS4</b>                                                                               |
| <b>Sensitivity</b>       | 5.6pg/mL                                                                                    |
| <b>Entrez</b>            | <b>3904</b>                                                                                 |
| <b>Expiration Date</b>   | Please enquire.                                                                             |