

Product Datasheet

Mouse LAMa1 ELISA Kit (orb440450)

Description	Mouse Laminin Alpha 1 (LAMa1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P19137
Sensitivity	0.29ng/mL
Entrez	16772
Expiration Date	Please enquire.