

Product Datasheet

Human LAT1 ELISA Kit (orb440466)

Description	Human Large Neutral Amino Acid Transporter 1 (LAT1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q01650
Sensitivity	0.056ng/mL
Entrez	8140
Expiration Date	Please enquire.