

Product Datasheet

Bovine LBP ELISA Kit (orb440470)

Description	Bovine Lipopolysaccharide Binding Protein (LBP) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q2TB10
Sensitivity	0.233ng/mL
Entrez	512242
Expiration Date	Please enquire.