

Product Datasheet

Human LDHB ELISA Kit (orb440485)

Description	Human Lactate Dehydrogenase B (LDHB) Sandwich ELISA Kit
Reactivity	Human
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P07195
Sensitivity	1.45ng/mL
Entrez	3945
Expiration Date	Please enquire.