

Product Datasheet

Human LECT2 ELISA Kit (orb440500)

Description	Human high sensitivity Leukocyte Cell Derived Chemotaxin 2 (LECT2) Sandwich ELISA Kit
Reactivity	Human
Range	0.23-15ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O14960
Sensitivity	0.05ng/mL
Entrez	3950
Expiration Date	Please enquire.