




Product Datasheet

CD47 antibody (orb44051)


Description	Mouse Monoclonal to CD47.
Reactivity	Human, Porcine, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, WB
Immunogen	COS-7 (African green monkey) cells
Target	CD47
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Non-reducing conditions.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Uniprot ID	Q08722
Entrez	961
Dilution Range	Flow cytometry: Recommended dilution: 2 µg/ml. Western blotting: Non-reducing conditions.
Expiration Date	12 months from date of receipt.



Flow cytometry (surface staining) of lys...



Separation of human lymphocytes (red-fil...



Flow cytometry surface staining pattern ...