

Product Datasheet

Chicken LEP ELISA Kit (orb440511)

Description	Chicken Leptin (LEP) Sandwich ELISA Kit
Reactivity	Gallus
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O42164
Sensitivity	14.8pg/mL
Entrez	373955
Expiration Date	Please enquire.