

Product Datasheet

Equine LEP ELISA Kit (orb440512)

Description	Equine Leptin (LEP) Sandwich ELISA Kit
Reactivity	Equine
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F7AZE4
Sensitivity	0.118ng/mL
Entrez	100034042
Expiration Date	Please enquire.