

Product Datasheet

Goat LEP ELISA Kit (orb440513)

Description	Goat Leptin (LEP) Sandwich ELISA Kit
Reactivity	Goat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q257X2
Sensitivity	0.081ng/mL
Entrez	100860794
Expiration Date	Please enquire.