

## Product Datasheet

### Human LEP ELISA Kit (orb440514)

|                          |                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Leptin (LEP) Sandwich ELISA Kit                                                                           |
| <b>Reactivity</b>        | Human                                                                                                           |
| <b>Range</b>             | 0.156-10ng/mL                                                                                                   |
| <b>Form/Appearance</b>   | Traditional                                                                                                     |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.                        |
| <b>Note</b>              | For research use only                                                                                           |
| <b>Application notes</b> | serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Sample Types</b>      | serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| <b>Uniprot ID</b>        | <b>P41159</b>                                                                                                   |
| <b>Sensitivity</b>       | 0.059ng/mL                                                                                                      |
| <b>Entrez</b>            | <b>3952</b>                                                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                                                                 |