

Product Datasheet

Rabbit LEP ELISA Kit (orb440518)

Description	Rabbit Leptin (LEP) Sandwich ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	B0KZL0
Sensitivity	0.054ng/mL
Entrez	100008747
Expiration Date	Please enquire.