

Product Datasheet

Rat LEPR ELISA Kit (orb440523)

Description	Rat Leptin Receptor (LEPR) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q62959
Sensitivity	0.66ng/mL
Entrez	24536
Expiration Date	Please enquire.