

Product Datasheet

Mouse LGALS3BP ELISA Kit (orb440525)

Description	Mouse Lectin Galactoside Binding, Soluble 3 Binding Protein (LGALS3BP) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q07797
Sensitivity	0.125ng/mL
Entrez	19039
Expiration Date	Please enquire.