

Product Datasheet

Bovine LIF ELISA Kit (orb440536)

Description	Bovine Leukemia Inhibitory Factor (LIF) Sandwich ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q27956
Sensitivity	11.5pg/mL
Entrez	280840
Expiration Date	Please enquire.